





Model 430 OEM Single Phase Lock-In Amplifier

Scitec Instruments' Model 430 OEM analogue lock-in amplifier uses advanced technology to provide a versatile high performance instrument.

- Single channel single phase operation
- Differential, single-ended or current input
- Gain settings from $3\mu\text{V}$ to 10V
- High performance wide bandwidth input gain stage
- Low pass filter output time constants from $100\mu\text{s}$ to 30s
- 1F and 2F reference signal operation
- 90° step and fine phase control

Input Modes: Differential, Single Ended, Current



INPUT SIGNAL CHANNELS

The input signal channel amplifies the input signal to a level suitable for the demodulator. High performance, low-noise, broad-band amplifiers are used throughout.

The input circuits can accept differential or single-ended inputs via the input BNC connectors. Jumper options within the unit allow the outer SMB contact or screen to act as a high impedance differential input or allow it to be connected to ground for single-ended operation. The 430 lock-in will also accept current inputs through the development of a voltage across a resistor. This mode is not guaranteed by Scitec.



The input channels are DC coupled rather than the more normal AC coupling seen on other lock-in amplifiers as the noise performance is improved.

Input	Differential or single ended voltage or current input via SMB socket
Sensitivity	3 μ V to 10V in 1, 3, 10 steps. The input gain is set using jumpers enabling simple gain changes. One set of gain resistors are mounted on solder pillars so that they can be easily changed by the user with the aid of a soldering iron.
Input Impedance	10 ¹² Ω 1nF, dc coupled
Frequency	10Hz to 100kHz
Maximum Inputs	\pm 10V before saturation occurs.
Noise	Scitec Instruments no longer specifies input noise values as this leads to comparison with other manufacturers data sheets which are clearly grossly in error. If you wish for details of these values then please contact us and we will explain the situation.
Gain Accuracy	1%
Gain Stability	200ppm/ $^{\circ}$ C
Dynamic Reserve	0dB to 80dB adjustable via jumpers.

DEMODULATOR

The input stage drives a high performance demodulator to recover the input signal.

OUTPUT

The demodulator output is passed through a low pass filter before being amplified for output.

Low Pass Filter Time Constant	100 μ S to 30s in 1, 3, 10 steps. The time constant is set using jumpers enabling simple time constant changes. One set of resistor capacitor values are mounted on solder pillars so that they can be easily modified with the aid of a soldering iron.
Outputs - SMB connectors	\pm 100mV to \pm 10V full scale output. Can be modified through jumper settings.

REFERENCE CHANNEL

The reference signal is used to generate the signals that drive the demodulator. A fine phase shifting circuit allows the reference signal to be phase shifted from 0 $^{\circ}$ to 150 $^{\circ}$ with relation to the signal input. A second circuit then produces signals that are phase shifted by 0 $^{\circ}$, 90 $^{\circ}$, 180 $^{\circ}$ and 270 $^{\circ}$ at both the reference frequency and twice the reference frequency. Any of these 8 signals can be used to drive the demodulator.

Frequency	10Hz to 100kHz
Trigger	Standard TTL, 95% mark/space ratio min. Rising edge triggered.
Acquisition time	10s max.
Phase control	90 $^{\circ}$ steps + fine shift in range 0 $^{\circ}$ - 150 $^{\circ}$
Phase Drift	0.1 $^{\circ}$ / $^{\circ}$ C
1F and 2F operation	

GENERAL

Power	-15V, 0V, +5V, +15V DC @ 35mA per supply. Power connections are made via 64 pin type C DIN41612 connector. Signal connections are not possible via this connector.
Mechanical	100 x 160mm
Temperature range	0-50 $^{\circ}$ C (operational)
Warranty	2 years from date of shipment

