








SPAD-8

8-Channel SPAD Module

- 8-channel SPAD detector module
- bh multi-dimensional TCSPC technique
- Interfaces directly to all bh TCSPC systems
- Simultaneous measurement in all 8 channels
- 1 x 8 arrangement of detector channels
- Instrument response width 70 ps FWHM
- Max. count rate > 5 MHz
- Thermo-electrically cooled
- Power supply and control via bh DCC-100 detector controller card



The SPAD-8 module contains eight actively quenched SPAD pixels on a single silicon chip. The signals of the SPADs are recorded by a single bh TCSPC module. The module uses bh's multi-dimensional TCSPC technique. For each photon, the SPAD-8 delivers a timing pulse and the number of the SPAD pixel that detected the photon. The TCSPC module builds up a photon distribution versus time and pixel number, or stores the individual events as time-tag data. The technique avoids any time gating or detector multiplexing and thus achieves a near-ideal counting efficiency. Power supply, SPAD excess-voltage control, an current for the TE cooler are provided by a bh DCC-100 detector controller card.



TE cooler current
SPAD excess voltage

Becker & Hickl GmbH
Nahmitzer Damm 30
12277 Berlin
Tel. +49 / 30 / 787 56 32
Fax. +49 / 30 / 787 57 34
email info@becker-hickl.com
http://www.becker-hickl.com



id Quantique
sales@idquantique.com
www@idquantique.com



Boston Electronics Corp
tcspc@boselec.com
www.boselec.com








SPAD-8

Specifications

Number of pixels	8
Pixel arrangement	1 by 8
Active area (each pixel)	40 x 40 μm
Pixel pitch, centre-centre	60 μm
Optical adapter	C mount
Spectral response	350 to 900 nm
Peak quantum efficiency, 500 nm	35 %
Channel uniformity	5%
Dark count rate, per channel	< 1000, TE cooler current 0.5A
IRF width, fwhm	70 ps (typical value)
Time skew between Channels	< 150 ps
Dead time	50 ns
Timing Output	SMA, 50 Ω , negative pulse
Routing Signal	3 bit + Error Signal, TTL/CMOS
Power Supply	From bh DCC-100 card
Dimensions	40 mm \times 40 mm \times 72 mm



Related Products

SPC-130 EM TCSPC modules
 SPC-150 TCSPC modules
 SPC-830 TCSPC modules
 SPC-630 TCSPC modules

Simple-Tau 130 compact TCSPC systems
 Simple-Tau 150 compact TCSPC systems
 Simple-Tau 830 compact TCSPC systems
 FLIM systems for laser scanning microscopy

DCC-100 detector controller
 PML-SPEC and MW-FLIM multi-wavelength detectors
 id-100 SPAD detector modules
 BDL-SMC and BHLF picosecond diode lasers

Related Literature

W. Becker, Advanced time-correlated single photon counting techniques. Springer 2005.
 W. Becker, The bh TCSPC Handbook, 466 pages, 503 references. Available on www.becker-hickl.com
 Please see also www.becker-hickl.com, 'Literature', 'Application notes'



More than 15 years experience in multi-dimensional TCSPC. More than 700 TCSPC systems worldwide.